

Afib

Double dose Gadavist

	Coverage	Gating	Breath Hold
Loc #1		none	NO
Loc #2: place heart at isocenter		none	NO
Twist Timer			
Axial ceMRA	Scan from skin to skin. Use single dose of contrast (dilute 50/50 with saline if you need more volume)and repeat if timing is wrong		end inspiration or end expiration
TI Scout 2 R-R	1 slice	prospective	YES
MDE single shot VLA	1 slice	prospective	optional
MDE single shot 4CH	1 slice	prospective	optional
MDE single shot LVOT	1 slice	prospective	optional
MDE single shot SAS	Left Ventricle	prospective	optional
TI Scout 1 R-R	1 slice	prospective	YES
3D MDE	from coronaries through ventricle	prospective	NO
IF 3D MDE is poor quality or ABNORMAL			
MDE segmented SAS	Left Ventricle	prospective	multiple
MDE segmented VLA	1 slice	prospective	YES
MDE segmented 4CH	1 slice	prospective	YES
MDE segmented LVOT	1 slice	prospective	YES